

Anti-Brca1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P38398

Item Number: CM0013583



Introduction

Rabbit polyclonal antibody against human BRCA1 for WB, IHC-P, IF/ICC, and ELISA. Recognizes human BRCA1; ~208 kDa. Suitable for DNA damage repair and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	BRCA1 Rabbit pAb
Remarks/Alias	IRIS; PSCP; BRCAI; BRCC1; FANCS; PNCA4; RNF53; BROVCA1; PPP1R53; BRCA1
Species	Human
GeneID (Human)	672
GeneID	672
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P38398
Protein Weight	208kDa
Shipping	Ice bag