

Anti-Brca2 Rabbit Polyclonal Antibody For Wb - P51587

Item Number: CM0019862



Introduction

High-quality rabbit polyclonal antibody targeting human BRCA2 (UniProt P51587). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for DNA damage repair, cancer biology, and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	BRCA2 Rabbit pAb
Remarks/Alias	FAD; FACD; FAD1; GLM3; BRCC2; FANCD; PNCA2; FANCD1; XRCC11; BROVCA2; BRCA2
Species	Human
Gene ID	675 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1600 of BRCA2 (NP_000050.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51587
Protein Weight	384kDa
Shipping	Ice bag