

Anti-Brcc2/Blid Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Izy5

Item Number: CM0026149

Introduction

Anti-BRCC2/BLID rabbit polyclonal antibody suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. For apoptosis and cell death studies.

[Learn More](#)

Parameter	Value
Product Name	BRCC2/BLID Rabbit pAb
Remarks/Alias	BRCC2; BRCC2/BLID
Species	Human
Gene ID (Human)	414899
Gene ID	414899
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of human BRCC2/BLID (NP_001001786.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZY5
Protein Weight	12kDa
Shipping	Ice bag