


Anti-Brd2 Monoclonal Antibody For Wb, Elisa - P25440

Item Number: CM0005322



Introduction

Rabbit monoclonal antibody targeting human BRD2 (bromodomain-containing protein 2), validated for Western blot and ELISA. Raised against residues 1-345, it detects the ~88 kDa chromatin reader involved in gene regulation and chromatin architecture.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] BRD2 Rabbit mAb
Remarks/Alias	FSH; NAT; RNF3; FSRG1; RING3; FSHRG1; D6S113E; O27.1.1; BRD2-IT1; BRD2
Species	Human
GeneID (human)	6046
GeneID	6046
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-345 of human BRD2 (NP_005095.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P25440
Protein Weight	88kDa
Shipping	Ice bag