


Anti-Brd2 Rabbit Polyclonal Antibody For Wb, Elisa - P25440

Item Number: CM0022829



Introduction

Rabbit polyclonal antibody against human BRD2, validated for Western blot and ELISA. Recognizes bromodomain-containing protein 2, a chromatin reader involved in transcriptional regulation and genome compartmentalization. Suitable for chromatin biology and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	BRD2 Rabbit pAb
Remarks/Alias	FSH; NAT; RNF3; FSRG1; RING3; FSHRG1; D6S113E; O27.1.1; BRD2-IT1; BRD2
Species	Human
Gene ID	6046
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human BRD2 (NP_005095.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P25440
Protein Weight	88kDa
Shipping	Ice bag