


Anti-Brd4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O60885

Item Number: CM0003776



Introduction

High-quality anti-BRD4 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Recognizes human BRD4, an acetyl-histone reader crucial for epigenetic memory and transcription regulation. Ideal for cancer and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	BRD4 Rabbit mAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; BRD4
Species	Human
GeneID (human)	23476
GeneID	23476
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human BRD4 (O60885).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O60885
Protein Weight	152kDa
Shipping	Ice bag