


Anti-Brd9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H8M2

Item Number: CM0028241



Introduction

Rabbit polyclonal antibody targeting human BRD9, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for chromatin remodeling and homologous recombination research.

[Learn More](#)

Parameter	Value
Product Name	BRD9 Rabbit pAb
Remarks/Alias	PRO9856; SMARCI2; LAVS3040; BRD9
Species	Human
GeneID (human)	65980
GeneID	65980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 274-597 of human BRD9 (NP_076413.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8M2
Protein Weight	67 kDa
Shipping	Ice bag