


Anti-Brf2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Haw0

Item Number: CM0026306



Introduction

Rabbit polyclonal antibody targeting human BRF2. Validated for WB and ELISA, recognizes a 47kDa band in human, mouse, and rat samples. Suitable for transcription and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	BRF2 Rabbit pAb
Remarks/Alias	BRFU; TFIIB50; BRF2
Species	Human
GeneID (human)	55290
GeneID	55290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-419 of human BRF2 (NP_060780.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HAW0
Protein Weight	47kDa
Shipping	Ice bag