

Anti-Brg1/Smarca4 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P51532

Item Number: CM0003899



Introduction

Monoclonal antibody targeting human BRG1/SMARCA4 (SWISS: P51532), validated for WB, IF/ICC, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for chromatin remodeling and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	BRG1/SMARCA4 Rabbit mAb
Remarks/Alias	BRG1; CSS4; SNF2; SWI2; MRD16; RTPS2; BAF190; SNF2L4; SNF2LB; hSNF2b; BAF190A; SNF2-beta; BRG1/SMARCA4
Species	Human
Gene ID (Human)	6597
Gene ID	6597
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human BRG1/SMARCA4 (P51532)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51532
Protein Weight	185 kDa
Shipping	Ice bag