

Anti-Brg1/Smarca4 Polyclonal Antibody For Wb, If/Icc, Elisa - P51532

Item Number: CM0014190



Introduction

Rabbit polyclonal antibody against human BRG1/SMARCA4 (P51532). Validated for western blot, immunofluorescence/immunocytochemistry, and ELISA. Chromatin remodeling factor involved in gene regulation and EMT.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] BRG1/SMARCA4 Rabbit pAb
Remarks/Alias	BRG1; CSS4; SNF2; SWI2; MRD16; RTP52; BAF190; SNF2L4; SNF2LB; hSNF2b; BAF190A; SNF2-beta; A4
Species	Human
Gene ID (Human)	6597
Gene ID	6597
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P51532
Protein Weight	185kDa
Shipping	Ice bag