

Anti-Brinp3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q76B58

Item Number: CM0021264

Introduction

Rabbit polyclonal antibody targeting BRINP3 (FAM5C). For WB, IF/ICC, ELISA applications. Reacts with human, mouse, rat. Recognizes 88 kDa protein implicated in neuronal proliferation inhibition and pituitary tumorigenesis. UniProt Q76B58.

[Learn More](#)

Parameter	Value
Product Name	BRINP3 Rabbit pAb
Remarks/Alias	FAM5C; DBCCR1L; DBCCR1L1; BRINP3
Species	Human
GeneID (human)	339479
GeneID	339479
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q76B58
Protein Weight	88kDa
Shipping	Ice bag