

Anti-Brip1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bx63

Item Number: CM0027587



Introduction

Rabbit polyclonal antibody against human BRIP1 (FANCJ). Validated for Western blot and ELISA. Useful for studying DNA repair, Fanconi anemia pathway, and chromosomal stability.

[Learn More](#)

Parameter	Value
Product Name	BRIP1 Rabbit pAb
Remarks/Alias	OF; BACH1; FANCJ; BRIP1
Species	Human
GeneID (Human)	83990
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human BRIP1 (NP_114432.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BX63
Protein Weight	141kDa
Shipping	Ice bag