


Anti-Brix1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Tdn6

Item Number: CM0013019



Introduction

Rabbit polyclonal antibody against human BRIX1, validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for ribosome biogenesis and nucleolar function studies.

[Learn More](#)

Parameter	Value
Product Name	BRIX1 Rabbit pAb
Remarks/Alias	BRIX; BXDC2; BRIX1
Species	Human
Gene ID (Human)	55299
Gene ID	55299
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human BRIX1 (NP_060791.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDN6
Protein Weight	41kDa
Shipping	Ice bag