

Anti-Brk1 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q8Wuw1

Item Number: CM0020634



Introduction

Rabbit polyclonal antibody targeting human BRK1 (Protein BRICK1, HSPC300). Validated for IHC-P, IF/ICC, and ELISA with cross-reactivity to human and mouse. Ideal for studying WAVE complex-mediated actin regulation and endocytic trafficking pathways.

[Learn More](#)

Parameter	Value
Product Name	BRK1 Rabbit pAb
Remarks/Alias	MDS027; hHBrk1; C3orf10; HSPC300; BRK1
Species	Human
Gene ID (Human)	55845
Gene ID	55845
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-75 of human BRK1 (NP_060932.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WUW1
Protein Weight	9kDa
Shipping	Ice bag