


Anti-Brms1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Hcu9

Item Number: CM0003271



Introduction

Recombinant rabbit monoclonal antibody targeting human BRMS1 for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for apoptosis, metastasis suppression, and NF-kB signaling research.

[Learn More](#)

Parameter	Value
Product Name	BRMS1 Rabbit mAb
Remarks/Alias	BRMS1; breast cancer metastasis-suppressor 1
Species	Human
Gene ID (Human)	25855
Gene ID	25855
Immunogen	Synthetic Peptide (a synthetic peptide corresponding to a sequence within amino acids 1-100 of human BRMS1 (Q9HCU9))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCU9
Protein Weight	28kDa
Shipping	Ice bag