


Anti-Brms1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Hcu9

Item Number: CM0015294



Introduction

High-quality rabbit polyclonal antibody against human BRMS1, a metastasis suppressor. Validated for WB, IHC-P, IF/ICC, ELISA. Suitable for breast carcinoma and anoikis research. Reacts with human, mouse, rat.

[Learn More](#)

Parameter	Specification
Product Name	BRMS1 Rabbit pAb
Remarks/Alias	BRMS1
Species	Human
Gene ID (Human)	25855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-210 of human BRMS1 (NP_056214.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCU9
Protein Weight	28kDa
Shipping	Ice bag