


Anti-Brpf1 Polyclonal Antibody For Wb, If/Icc, Elisa - P55201

Item Number: CM0022355



Introduction

BRPF1 Rabbit Polyclonal Antibody (Peregrin) validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	BRPF1 Rabbit pAb
Remarks/Alias	BR140; IDDDFP; BRPF1
Species	Human
Gene ID (Human)	7862
Gene ID	7862
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human BRPF1 (NP_001003694)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55201
Protein Weight	137kDa
Shipping	Ice bag