

Anti-Brrf2 Monoclonal Antibody For Wb, If/Icc, Elisa - P03210

Item Number: CM0007673

Introduction

Rabbit monoclonal antibody targeting Epstein-Barr virus tegument protein BRRF2. Suitable for WB, IF/ICC, and ELISA. Ideal for studying viral immune evasion and host-pathogen interactions. Uses recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	BRRF2 Rabbit mAb
Remarks/Alias	Tegument protein BRRF2; BRRF2
Species	Epstein-Barr virus (strain B95-8) (HHV-4)
Gene ID	3783729
Immunogen	Recombinant protein of Epstein-Barr virus (strain B95-8) (HHV-4) BRRF2.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Epstein-Barr virus
SWISS	P03210
Protein Weight	56KDa
Shipping	Ice bag