

Anti-Brrf2 Rabbit Polyclonal Antibody For Wb, Elisa - P03210

Item Number: CM0024054

Introduction

Rabbit polyclonal antibody against Epstein-Barr virus tegument protein BRRF2, validated for WB and ELISA with human cross-reactivity. Ideal for studying viral immune evasion and cGAS-DNA phase separation inhibition.

[Learn More](#)

Parameter	Value
Product Name	BRRF2 Rabbit pAb
Remarks/Alias	Tegument protein BRRF2; BRRF2
Species	Epstein-Barr virus (strain B95-8) (HHV-4)
GeneID	3783729
Immunogen	Recombinant protein of Epstein-Barr virus (strain B95-8) (HHV-4) BRRF2
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P03210
Protein Weight	56kDa
Shipping	Ice bag