

Anti-Brsk1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Tdc3

Item Number: CM0024692

Introduction

Anti-BRSK1 rabbit polyclonal antibody for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. Ideal for studying neuronal polarization, centrosome duplication, and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	BRSK1 Rabbit pAb
Remarks/Alias	SAD-B; hSAD1; BRSK1
Species	Human
Gene ID (Human)	84446
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-510 of human BRSK1 (NP_115806.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDC3
Protein Weight	85 kDa
Shipping	Ice bag