


Anti-Brwd1 Polyclonal Antibody For Wb, Elisa - Q9Nsi6

Item Number: CM0013203



Introduction

Rabbit polyclonal antibody targeting BRWD1 (WDR9, DCAF19), suitable for WB and ELISA. Cross-reacts with human and mouse proteins. Ideal for investigating chromatin remodeling, transcriptional activation, and ciliary biology.

[Learn More](#)

Parameter	Value
Product Name	BRWD1 Rabbit pAb
Remarks/Alias	N143; WDR9; WRD9; DCAF19; C21orf107; BRWD1
Species	Human
GeneID (Human)	54014
GeneID	54014
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human BRWD1 (NP_001007247.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NSI6
Protein Weight	263kDa
Shipping	Ice bag