

Anti-Bsn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Upa5

Item Number: CM0022858

Introduction

High-quality BSN rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects endogenous bassoon in human and mouse brain tissue. Essential for synaptic biology and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	BSN Rabbit pAb
Remarks/Alias	ZNF231; BSN
Species	Human
Gene ID	8927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 830-950 of human BSN (NP_003449.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UPA5
Protein Weight	416kDa
Shipping	Ice bag