

Anti-Bsnd Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wz55

Item Number: CM0024862

Introduction

Rabbit polyclonal antibody targeting BSND (Barttin), a chloride channel regulator, for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for renal and auditory research.

[Learn More](#)

Parameter	Specification
Product Name	BSND Rabbit pAb
Remarks/Alias	BART; DFNB73; BSND
Species	Human
Gene ID (Human)	7809
Gene ID	7809
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 54-320 of human BSND (NP_476517.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WZ55
Protein Weight	35kDa
Shipping	Ice bag