

Anti-Bst2 (Cd317/Tetherin) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q10589

Item Number: CM0017709



Introduction

Rabbit polyclonal antibody against human BST2 (CD317/Tetherin) validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human and mouse BST2. Ideal for antiviral restriction, immune signaling, and multiple myeloma studies.

[Learn More](#)

Parameter	Value
Product Name	BST2 Rabbit pAb
Remarks/Alias	CD317; HM1.24; TETHERIN; BST2
Species	Human
GeneID (Human)	684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-161 of human BST2 (NP_004326.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q10589
Protein Weight	20kDa
Shipping	Ice bag