

Anti-Bst2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q10589

Item Number: CM0007383



Introduction

CamelBio offers an Anti-BST2 Rabbit Monoclonal Antibody (Q10589) targeting human Tetherin/CD317. Suitable for WB, IHC-P, IF/ICC, and ELISA. Ideal for studying viral restriction and multiple myeloma.

[Learn More](#)

Parameter	Value
Product Name	BST2 Rabbit mAb
Remarks/Alias	CD317; HM1.24; TETHERIN; BST2
Species	Human
Gene ID (Human)	684
Gene ID	684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-160 of human BST2 (NP_004326.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q10589
Protein Weight	20kDa
Shipping	Ice bag