

Anti-Bsu1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Lr78

Item Number: CM0023021

Introduction

Rabbit polyclonal antibody targeting Arabidopsis thaliana BSU1, a positive regulator of brassinosteroid signaling. Validated for Western blot (WB) and ELISA. Ideal for plant hormone signaling and growth development research.

[Learn More](#)

Parameter	Value
Product Name	BSU1 Rabbit pAb
Remarks/Alias	BRI1 SUPPRESSOR 1; F21B7.7; BSU1
Species	Arabidopsis thaliana
Gene ID	838804
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 694-793 of Arabidopsis thaliana BSU1 (NP_171844.6)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q9LR78
Protein Weight	88kDa
Shipping	Ice bag