

Anti-Btbd7 Rabbit Polyclonal Antibody For Wb, Elisa - Q9P203

Item Number: CM0023828

Introduction

Anti-BTBD7 Rabbit Polyclonal Antibody validated for Western blot and ELISA, detecting BTBD7 in Human, Mouse, and Rat samples. A key tool for epithelial dynamics and organ branching research in IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	BTBD7 Rabbit pAb
Remarks/Alias	FUP1; BTBD7
Species	Human
Gene ID (Human)	55727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 5-75 of human BTBD7 (NP_060637.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P203
Protein Weight	126kDa
Shipping	Ice bag