

Anti-Btbd9 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Q07

Item Number: CM0019829

Introduction

The BTBD9 Rabbit pAb, raised against amino acids 283-582 of human BTBD9, detects endogenous BTBD9 in Western blot and ELISA. Ideal for neuroscience and kidney research.

[Learn More](#)

Parameter	Value
Product Name	BTBD9 Rabbit pAb
Remarks/Alias	dj322i12.1; BTBD9
Species	Human
GeneID (Human)	114781
GeneID	114781
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 283-582 of human BTBD9 (AAI01358.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96Q07
Protein Weight	69kDa
Shipping	Ice bag