

Anti-Btc Rabbit Polyclonal Antibody For Wb, Elisa - P35070

Item Number: CM0023013



Introduction

Polyclonal antibody targeting human BTC (Probetacellulin). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for growth factor signaling studies in cancer and development research.

[Learn More](#)

Parameter	Value
Product Name	BTC Rabbit pAb
Remarks/Alias	BTC
Species	Human
GeneID (Human)	685
GeneID	685
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35070
Protein Weight	20kDa
Shipping	Ice bag