


Anti-Btf3 Rabbit Polyclonal Antibody For Wb, Elisa - P20290

Item Number: CM0023545



Introduction

BTF3 rabbit polyclonal antibody validated for Western blot and ELISA in human, mouse, and rat. BTF3 (NAC-beta) functions as a nascent polypeptide-associated complex subunit and transcription factor.

[Learn More](#)

Parameter	Value
Product Name	BTF3 Rabbit pAb
Remarks/Alias	NACB; BTF3a; BTF3b; BETA-NAC; BTF3
Species	Human
GeneID (Human)	689
GeneID	689
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 50-150 of human BTF3 (NP_001198.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20290
Protein Weight	22kDa
Shipping	Ice bag