

Anti-Btg1 Polyclonal Antibody For Wb, If/Icc, Elisa - P62324

Item Number: CM0025223

Introduction

BTG1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA. Reacts with human BTG1, an anti-proliferative protein implicated in cancer. Suitable for research on cell cycle regulation and chromosomal rearrangements.

[Learn More](#)

Parameter	Value
Product Name	BTG1 Rabbit pAb
Remarks/Alias	APRO2; BTG1
Species	Human
Gene ID (Human)	694
Gene ID	694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-171 of human BTG1 (NP_001722.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P62324
Protein Weight	19kDa
Shipping	Ice bag