

Anti-Btg3 Polyclonal Antibody For Wb, Elisa - Q14201

Item Number: CM0025521

Introduction

Anti-BTG3 Rabbit Polyclonal Antibody, validated in WB and ELISA, targets human BTG3 (Q14201), a cell cycle regulator abundant in neuroepithelium and reproductive tissues. Ideal for cancer and developmental research.

[Learn More](#)

Parameter	Value
Product Name	BTG3 Rabbit pAb
Remarks/Alias	ANA; TOB5; TOFA; APRO4; TOB55; ANA/BTG3; BTG3
Species	Human
Gene ID (Human)	10950
Gene ID	10950
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 150-250 of human BTG3 (NP_006797.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14201
Protein Weight	29kDa
Shipping	Ice bag