


Anti-Btk Polyclonal Antibody For Wb, Ihc-P, Elisa - Q06187

Item Number: CM0018107



Introduction

BTK rabbit polyclonal antibody suited for Western blot, IHC-P, and ELISA applications. Recognizes human, mouse, and rat BTK, a critical kinase in B-cell development and signaling. Product referenced to UniProt Q06187.

[Learn More](#)

Parameter	Value
Product Name	BTK Rabbit pAb
Remarks/Alias	AT; ATK; BPK; XLA; IMD1; AGMX1; IGHD3; PSCTK1; BTK
Species	Human
GeneID	695
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-240 of human BTK (NP_000052.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06187
Protein Weight	76kDa
Shipping	Ice bag