


Anti-Btk Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q06187

Item Number: CM0002078



Introduction

BTK Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, ELISA with cross-reactivity to human and mouse. Ideal for studying B-cell receptor signaling and immunodeficiency research.

[Learn More](#)

Parameter	Value
Product Name	BTK Rabbit mAb
Remarks/Alias	AT; ATK; BPK; XLA; IMD1; AGMX1; IGHD3; PSCTK1; BTK
Species	Human
GeneID (Human)	695
GeneID	695
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 560-659 of human BTK (Q06187)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06187
Protein Weight	76 kDa
Shipping	Ice bag