

Anti-Btla (Cd272) Polyclonal Antibody For Wb, Elisa - Q7Z6A9

Item Number: CM0017589



Introduction

Rabbit polyclonal antibody against human B- and T-lymphocyte attenuator (BTLA, CD272) for Western blot (WB) and ELISA. Detects human and mouse BTLA. Suitable for immune checkpoint and lymphocyte signaling studies.

[Learn More](#)

Parameter	Value
Product Name	BTLA Rabbit pAb
Remarks/Alias	BTLA1; CD272; BTLA
Species	Human
Gene ID (Human)	151888
Gene ID	151888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-241 of human BTLA (NP_001078826.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z6A9
Protein Weight	33kDa
Shipping	Ice bag