

Anti-Btla Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z6A9

Item Number: CM0027408



Introduction

BTLA (CD272) rabbit polyclonal antibody, validated for Western blot and ELISA, cross-reactive with mouse and rat. Ideal for studying lymphocyte attenuation and immune checkpoint signaling. Protein weight: 33 kDa.

[Learn More](#)

Parameter	Value
Product Name	BTLA Rabbit pAb
Remarks/Alias	BTLA1, CD272, BTLA
Species	Human
GeneID (Human)	151888
GeneID	151888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-241 of human BTLA (NP_001078826.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q7Z6A9
Protein Weight	33kDa
Shipping	Ice bag