

Anti-Btn2A1 Rabbit Monoclonal Antibody For Wb, Elisa - Q7Kyr7

Item Number: CM0007340



Introduction

Recombinant rabbit monoclonal antibody targeting human butyrophilin subfamily 2 member A1 (BTN2A1). Suitable for Western blot and ELISA; recognizes native and recombinant human BTN2A1. Ideal for immunology and membrane protein studies.

[Learn More](#)

Parameter	Value
Product Name	BTN2A1 Rabbit mAb
Remarks/Alias	BTF1; BT2.1; BTN2.1; DJ3E1.1; BK14H9.1; BTN2A1
Species	Human
GeneID (human)	11120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human BTN2A1 (NP_008980.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7KYR7
Protein Weight	60kDa
Shipping	Ice bag