

Anti-Btn2A1 Rabbit Polyclonal Antibody For Wb And Elisa - Q7Kyr7

Item Number: CM0016404



Introduction

Rabbit polyclonal antibody targeting human butyrophilin subfamily 2 member A1 (BTN2A1). Suitable for Western blot and ELISA. Ideal for cell surface receptor and tissue expression studies.

[Learn More](#)

Parameter	Value
Product Name	BTN2A1 Rabbit pAb
Remarks/Alias	BTF1; BT2.1; BTN2.1; DJ3E1.1; BK14H9.1; BTN2A1
Species	Human
Gene ID	11120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-240 of human BTN2A1 (NP_510961.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7KYR7
Protein Weight	60kDa
Shipping	Ice bag