

Anti-Btn3A1/2/3 Monoclonal Antibody For Fc - 000481 / P78410 / 000478

Item Number: CM0010074



Introduction

Rabbit anti-Human BTN3A1/2/3 mAb for flow cytometry. Targets butyrophilin subfamily 3 members involved in T-cell activation and immune response. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human BTN3A1/2/3 mAb
Remarks/Alias	BTf5; BT3.1; CD277; BTN3.1
Species	Human
GeneID (Human)	11119/11118/10384
GeneID	11119/11118/10384
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 470-630 of human BTN3A1 (NP_008979.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O00481/P78410/O00478
Protein Weight	58kDa
Shipping	Ice bag