

Anti-Btn3A3 Rabbit Polyclonal Antibody For Wb, Elisa - O00478

Item Number: CM0018063



Introduction

BTN3A3 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat BTN3A3, a T-cell adaptive immune response regulator. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	BTN3A3 Rabbit pAb
Remarks/Alias	BTF3; BTN3.3; BTN3A3
Species	Human
GeneID (Human)	10384
GeneID	10384
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-584 of human BTN3A3 (NP_008925.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00478
Protein Weight	65kDa
Shipping	Ice bag