

Anti-Btnl2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uir0

Item Number: CM0025673

Introduction

BTNL2 rabbit polyclonal antibody (pAb) validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Detects human, mouse, and rat BTNL2, a negative regulator of T-cell proliferation. Suitable for immunology and immuno-oncology research.

[Learn More](#)

Parameter	Specification
Product Name	BTNL2 Rabbit pAb
Remarks/Alias	SS2, BTN7, BTL-II, HSBLMHC1, BTNL2
Species	Human
GeneID	56244
Immunogen	Recombinant fusion protein containing a.a. 246-455 of human BTNL2 (NP_062548.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIR0
Protein Weight	50kDa
Shipping	Ice bag