


Anti-Bub1 Rabbit Polyclonal Antibody For Wb, Elisa - 043683

Item Number: CM0024006



Introduction

Rabbit polyclonal antibody targeting human BUB1, validated for WB and ELISA. Detects endogenous BUB1 in human and mouse samples. Ideal for spindle assembly checkpoint and chromosome alignment research.

[Learn More](#)

Parameter	Value
Product Name	BUB1 Rabbit pAb
Remarks/Alias	BUB1A; BUB1L; hBUB1; MCPH30; BUB1
Species	Human
GeneID (Human)	699
GeneID	699
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human BUB1 (NP_004327.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43683
Protein Weight	122kDa
Shipping	Ice bag