

Anti-Bub1B Polyclonal Antibody For Wb, If/Icc, Elisa - O60566

Item Number: CM0023837



Introduction

Rabbit polyclonal antibody against BUB1B (BUBR1), validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for mitotic checkpoint and spindle assembly checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	BUB1B Rabbit pAb
Remarks/Alias	MVA1; SSK1; BUBR1; Bub1A; MAD3L; hBUBR1; BUB1beta; BUB1B
Species	Human
GeneID	701
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human BUB1B (NP_001202.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60566
Protein Weight	120kDa
Shipping	Ice bag