

Anti-Bub3 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43684

Item Number: CM0011120

Introduction

High-quality rabbit monoclonal antibody against human BUB3 (O43684). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for mitotic checkpoint and kinetochore research.

[Learn More](#)

Parameter	Value
Product Name	BUB3 Rabbit mAb
Remarks/Alias	BUB3L; hBUB3; BUB3
Species	Human
Gene ID (Human)	9184
Gene ID	9184
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human BUB3 (O43684).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43684
Protein Weight	37kDa
Shipping	Ice bag