

Anti-Bub3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43684

Item Number: CM0027139

Introduction

Rabbit polyclonal antibody targeting human BUB3 (37 kDa) mitotic checkpoint protein. Validated in WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. For spindle assembly checkpoint and chromosome segregation research.

[Learn More](#)

Parameter	Value
Product Name	BUB3 Rabbit pAb
Remarks/Alias	BUB3L; hBUB3; BUB3
Species	Human
GeneID	9184
Immunogen	Recombinant fusion protein containing amino acids 1-328 of human BUB3 (NP_004716.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43684
Protein Weight	37kDa
Shipping	Ice bag