


## Anti-Butyrylcholinesterase Monoclonal Antibody For Wb And Elisa - P06276

Item Number: CM0011178



### Introduction

This rabbit monoclonal antibody targets human Butyrylcholinesterase (P06276). Validated for WB and ELISA, it is suitable for research on neurotransmitter degradation and organophosphate detoxification.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Butyrylcholinesterase Rabbit mAb                                                                                                       |
| Remarks/Alias    | E1; CHE1; CHE2; BCHED; Butyrylcholinesterase                                                                                           |
| Species          | Human                                                                                                                                  |
| GeneID (Human)   | 590                                                                                                                                    |
| GeneID           | 590                                                                                                                                    |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 500-602 of human Butyrylcholinesterase (P06276). |
| Source           | Rabbit                                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                                  |
| Application      | WB, ELISA                                                                                                                              |
| Cross Reactivity | Human                                                                                                                                  |
| SWISS            | P06276                                                                                                                                 |
| Protein Weight   | 68kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |