


Anti-Butyrylcholinesterase Rabbit Polyclonal Antibody For Wb, Elisa - P06276

Item Number: CM0018494



Introduction

Rabbit polyclonal antibody targeting human Butyrylcholinesterase (BCHE). Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Supports neurobiology and toxicology research and assay development.

[Learn More](#)

Parameter	Value
Product Name	Butyrylcholinesterase Rabbit pAb
Remarks/Alias	E1; CHE1; CHE2; BCHED; Butyrylcholinesterase
Species	Human
Gene ID (Human)	590
Gene ID	590
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06276
Protein Weight	68kDa
Shipping	Ice bag