

Anti-Bves (Popdc1) Polyclonal Antibody For Wb And Elisa - Q8Ne79

Item Number: CM0015504

Introduction

Rabbit polyclonal antibody targeting human BVES (POPDC1), validated for Western blot and ELISA. Suitable for research on tight junctions, cardiac function, and muscle regeneration. Cross-reactivity: human. Predicted MW: 41 kDa.

[Learn More](#)

Parameter	Value
Product Name	BVES Rabbit pAb
Remarks/Alias	POP1; HBVES; CARICK; LGMD2X; POPDC1; LGMDR25; BVES
Species	Human
Gene ID	11149
Immunogen	Recombinant fusion protein containing amino acids 241-360 of human BVES (NP_001186492.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NE79
Protein Weight	41kDa
Shipping	Ice bag