


## Anti-Bysl Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13895

Item Number: CM0016842



### Introduction

Validated BYSL rabbit polyclonal antibody for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets bystin, a 50 kDa protein involved in rRNA processing and ribosome biogenesis. Ideal for placental and adhesion studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | BYSL Rabbit pAb                                                                                                   |
| Remarks/Alias    | Enp1; BYSTIN; BYSL                                                                                                |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 705                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 120-280 of human BYSL (NP_004044.3) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q13895                                                                                                            |
| Protein Weight   | 50kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |