

Anti-Bzip46 Polyclonal Antibody For Wb, Elisa - Q69Tw5

Item Number: CM0023957

Introduction

Rabbit polyclonal antibody against rice bZIP46 transcription factor, validated for WB and ELISA. Suitable for abscisic acid signaling and stress response research.

[Learn More](#)

Parameter	Value
Product Name	bZIP46 Rabbit pAb
Remarks/Alias	P0701E03.2; bZIP46
Species	Oryza sativa
GeneID	4340462
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of Oryza sativa bZIP46. (Q69TW5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q69TW5
Protein Weight	35kDa
Shipping	Ice bag